



Group 1 & 2 Seating Guide

REV 1 June 2026



SLING	STADIUM / FISH ON	CAPTAIN SEAT	PAN SEAT- CAPTAIN BACKREST	REHAB FRAME PAN SEAT WITH TA BACKREST
 DESCRIPTION Flexible fabric seat and backrest KEY FEATURES Foldable, great for traveling and includes 7/8" diameter back canes COMPATIBLE CHAIR MODELS Travel-Ease series Standard & Heavy Duty (P101, P181, P182, P183) SEAT WIDTH Standard sizes: 16", 18", 20" HD sizes: 22", 24", 26" SEAT DEPTH Standard sizes: 16" *non-adjustable HD sizes: 18", 20" *non-adjustable BACKREST HEIGHT Standard sizes: 16", 18", 20" HD sizes: 22", 24", 26"	 DESCRIPTION Contoured cushion seat with a short contoured backrest KEY FEATURES Great combo of sleek + comfort COMPATIBLE CHAIR MODELS EZ-GO (P321B) Roadster 3 Roadster 4 Roadster S3 Roadster S4 SEAT WIDTH 18" SEAT DEPTH 16" *non-adjustable BACKREST HEIGHT 18"	 DESCRIPTION Contoured cushion seat with a tall contoured backrest KEY FEATURES Ultimate comfort, full back support, lumbar & lateral contours COMPATIBLE CHAIR MODELS Vision CF (P322A) Vision Sport no lift or with lift (P326E, P326D) Vision Ultra (P325) Regal (P310) Dualer (P312) Gemini no lift or with lift (P301, P3011) Vision Super no lift or with lift (P327, P3274) Atlantis (P710) SEAT WIDTH Standard sizes: 16", 18", 20" HD sizes: 20", 22", 24", 26" SEAT DEPTH RANGE Standard sizes: 16" - 20" HD sizes: 18" - 22" *adjustable in 1" increments BACKREST HEIGHT Standard sizes: 18", 20" HD sizes: 20", 22"	 DESCRIPTION Solid seat base with a captain-style backrest KEY FEATURES Customizable sizing for width and depth as well as cushion thickness COMPATIBLE CHAIR MODELS Vision CF (P322A) Vision Sport no lift or with lift (P326E, P326D) Vision Ultra (P325) Regal (P3101) Dualer (P312) Gemini no lift or with lift (P301, P30111) Vision Super no lift or with lift (P327, P3274) Atlantis (P710) SEAT WIDTH Standard sizes: 16", 18", 20" HD sizes: 20", 22", 24", 26" SEAT DEPTH RANGE Standard sizes: 16" - 20" HD sizes: 18" - 22" *adjustable in 1" increments BACKREST HEIGHT Standard sizes: 18", 20" HD sizes: 20", 22" Usually ordered with a cushion (See PC Cushion Line)	 DESCRIPTION Rehab frame seat base with tension-adjustable backrest straps. R600 for one-seat post units R601 for four-seat posts units KEY FEATURES Adjustable backrest and multiple seat cushion options to support various spine & hip conditions. Includes 7/8" diameter back canes. COMPATIBLE CHAIR MODELS R600: Vision Ultra (P325) Vision Sport no lift (P326E) R601: Vision Super no lift (P327) SEAT WIDTH R600 sizes: 16", 18", 20", 22" R601 sizes: 18", 20", 22", 24" SEAT DEPTH RANGE R600 sizes: 16" - 20" R601 sizes: 18" - 22" *adjustable in 1" increments BACKREST HEIGHT R600 sizes: 18" R601 sizes: 18" Usually ordered with a cushion (See PC Cushion Line)

All specifications and information contained herein are believed to be accurate at the time of publication and are subject to change without notice. Not all sizes available on all models, please refer to each wheelchair's order form for available size options.

Enhancing Life Through Mobility.

©Merits Health Products, Inc. | 4245 Evans Avenue Fort Myers, FL 33901 | www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



Group 2 & 3 Seating Guide

REV 1 June 2026



CONTOURED CAPTAIN

DESCRIPTION

Entry-level seating with a contoured cushion seat and a tall contoured backrest.

Ultimate comfort, full back support, lumbar & lateral contours.

KEY FEATURES

- Compatible with VR2 electronics
- Compatible with captain-style headrest
- Compatible with elevate power function only

COMPATIBLE CHAIR MODELS

Ultra (P325)
Ultra CST (P325I)
Velocity (P325)
Vector HD (P323)

SEAT WIDTH

Standard sizes: 16", 18", 20", 22"
HD sizes: 20", 22", 24"

SEAT DEPTH

*adjustable in 1" increments
Standard sizes: 16" - 22"
HD sizes: 20" - 22"

BACKREST HEIGHT

Standard sizes: 18", 20", 22"
HD sizes: 20", 22"

PAN SEAT WITH CAPTAIN BACKREST

DESCRIPTION

Entry-level seating with a solid seat pan and a captain-style backrest.

Allows for customizable seat cushion width and depth.

KEY FEATURES

- Compatible with VR2 electronics
- Compatible with captain-style headrest
- Compatible with elevate power function only

COMPATIBLE CHAIR MODELS

Ultra (P325)
Velocity (P325)
Vector HD (P323)

SEAT WIDTH

Standard sizes: 16", 18", 20"
HD sizes: 22", 24"

SEAT DEPTH

*adjustable in 1" increments
Standard sizes: 16" - 20"
HD sizes: 18" - 22"

BACKREST HEIGHT

Standard sizes: 18", 20"
HD sizes: 20", 22"

Usually ordered with a cushion
(See PC Cushion Line)
Commonly ordered with the
PC Lite cushion

POWERED REHAB FRAME PAN SEAT AVID SEATING SYSTEM

DESCRIPTION

Mid-level power seating system with a solid pan seat and a solid pan backrest.

Allows for customizable seat cushion width and depth.

KEY FEATURES

- Compatible with VR2 or R-net electronics
- Compatible with midline or multi-axis headrest
- Compatible with trunk & pelvic supports
- Compatible with all power functions: tilt, recline, elevate, and leg rest options
- *No recline on Ultra models
- Standardly comes with back canes

COMPATIBLE CHAIR MODELS

Ultra (P325)
Ultra X (P325X)
Ultra 3A (P325)
Ultra HD (P331)
Velocity (P325)
Vector HD (P323)

SEAT WIDTH

Standard sizes: 16", 17", 18", 19", 20", 21", 22"
HD sizes: 18", 19", 20", 21", 22", 23", 24"

SEAT DEPTH RANGE

Standard sizes: 16", 17", 18", 19", 20", 21", 22"
HD sizes: 18", 19", 20", 21", 22", 23", 24"

BACKREST HEIGHT

Full & HD sizes: 17", 18", 22"
• Non-coded backrest: 17.5"
• Coded backrest: 18", 22"

Usually ordered with a cushion

POWERED REHAB FRAME PAN SEAT FREQUENCY SEATING SYSTEM

DESCRIPTION

Full-featured power seating system with a solid pan seat and a solid pan backrest with shear reduction.

Allows for customizable seat cushion width and depth.

KEY FEATURES

- 6" recline shear reduction
- Compatible with R-net electronics
- Compatible with midline or multi-axis headrest
- Compatible with trunk & pelvic supports
- Compatible with all power functions (tilt, recline, elevate, and leg rest options)

COMPATIBLE CHAIR MODELS

Axcel (P335)
Axcel Narrow (P335N)

SEAT WIDTH

Standard sizes: 14", 15", 16", 17", 18", 19", 20", 21", 22"

SEAT DEPTH RANGE

Standard sizes: 14", 15", 16", 17", 18", 19", 20", 21", 22"

BACKREST HEIGHT

Full & HD sizes: 17", 18", 22"
• Non-coded backrest: 17.5"
• Coded backrest: 18", 22"

Usually ordered with a cushion
(See PC Cushion Line)

*Not all sizes available on all models, please refer to each wheelchair's order form for available size options.

*All specifications and information contained herein are believed to be accurate at the time of publication and are subject to change without notice.

Enhancing Life Through Mobility.

©Merits Health Products, Inc. | 4245 Evans Avenue Fort Myers, FL 33901 | www.meritsusa.com | 239-772-0579 | sales@meritsusa.com