

merits

# ULTRA HD P331

Heavy Duty MWD Powerbase

K0837



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features

- Mid-wheel drive design provides outstanding maneuverability
- Adjustable seat width 20"-24" and seat depth 16"-20" for custom fit
- Optimized seating stability with 4-post seat design
- Tool-free battery box provides easy maintenance
- Front and rear suspension offers all day comfort and easy obstacle clearance
- Seat tilt angle-adjustable back 90° - 110°
- Dual in-line motors for enhanced efficiency, torque, range, and performance
- Unoccupied transit tie down brackets included

## Accessories

- USB charger • Oxygen tank holder bag

## Options

- Tilt 0° - 50°
- Articulating Footplate
- VR2 electronics
- Swing away mounting brackets and mounts
- Multiple Joystick options
- Power Seat Tilt/  
Angle-Adjustable Back
- Seat to Floor Height options
- Positioning Components  
and Headrests

## Specifications

|                           |                                                                                                                                                                            |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product#                  | P331                                                                                                                                                                       |
| Weight Capacity           | 450 lbs                                                                                                                                                                    |
| Drivetrain                | MWD                                                                                                                                                                        |
| Maximum Speed             | 5 mph                                                                                                                                                                      |
| Per-Charge Range          | up to 8 miles / 19 miles                                                                                                                                                   |
| Battery Size (2 required) | 55Ah, 22NF AGM<br>51Ah, 22NF GEL<br>33Ah, U1 AGM                                                                                                                           |
| Ground Clearance          | 2.875"                                                                                                                                                                     |
| Turning Radius            | 21"                                                                                                                                                                        |
| Overall Length            | 35" / 45" with footplate                                                                                                                                                   |
| Overall Width             | 27"                                                                                                                                                                        |
| Overall Height            | 42"                                                                                                                                                                        |
| Seat-to-Ground Height     | 17.5" - 18.5"                                                                                                                                                              |
| Seat Width                | 20" - 24"                                                                                                                                                                  |
| Seat Depth                | 16" - 20"                                                                                                                                                                  |
| Backrest Height           | 21" - 23"                                                                                                                                                                  |
| Total Weight              | 288lbs. - 310 lbs                                                                                                                                                          |
| Motor                     | 2-Pole/200W                                                                                                                                                                |
| Controller                | PG VR2-A2 60A                                                                                                                                                              |
| Charger                   | 6A off board                                                                                                                                                               |
| Gradient                  | 6°                                                                                                                                                                         |
| Front Wheels              | 6" x 2"                                                                                                                                                                    |
| Drive Wheels              | 12" x 3"                                                                                                                                                                   |
| Rear Wheels               | 6" x 2"                                                                                                                                                                    |
| HCPCSCode                 | K0837                                                                                                                                                                      |
| Color Options             | <input type="checkbox"/> White<br><input checked="" type="checkbox"/> Red<br><input checked="" type="checkbox"/> Blue<br><input checked="" type="checkbox"/> Charcoal Gray |



Features MK batteries



merits



**merits** Enhancing life through mobility.

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0825