

merits

# VISION SUPER with LIFT P3274

Super Heavy Duty MWD Powerbase

K0824 / K0825 / K0830 / K0831



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features

- True mid-wheel: 21" turning radius
- Elevating 10" power scissor seat
- Dynamic LiNX & Shark controller for convenient Bluetooth programming and diagnostics
- Premium, semi-reclining ergonomic seat up to 135°
- 18" - 22" customizable seat depth

## Accessories

- Unoccupied transit tie down kit
- Cup holder
- USB charger
- Oxygen tank holder
- Walker holder
- Foldable rear basket
- Safety flag
- Side under, arm bag

## Options

- Optional, 20" - 24" seat width
- Elevating leg rests (K0195)
- Captain or Pan powered seat option



## Specifications

|                           |                                                                                                                                                                                          |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product#                  | P3274 w/lift                                                                                                                                                                             |
| Weight Capacity           | 400lbs. w/lift                                                                                                                                                                           |
| Drivetrain                | MWD                                                                                                                                                                                      |
| Seat Type                 | Captain or Pan                                                                                                                                                                           |
| Seat Width                | 20" - 24"                                                                                                                                                                                |
| Seat Depth                | 18" - 22"                                                                                                                                                                                |
| Backrest Height           | 20" - 22"                                                                                                                                                                                |
| Seat-to-Floor Height      | 24" - 34" w lift                                                                                                                                                                         |
| Seat Elevation            | 10" - P3274 w/lift                                                                                                                                                                       |
| Overall Length            | 35.4"                                                                                                                                                                                    |
| Overall Width             | 27.2"                                                                                                                                                                                    |
| Overall Height            | 48"                                                                                                                                                                                      |
| Turning Radius            | 21"                                                                                                                                                                                      |
| Ground Clearance          | 2.2"                                                                                                                                                                                     |
| Total Weight              | 243 lbs.                                                                                                                                                                                 |
| Wt of Heaviest Pc.        | 121 lbs.                                                                                                                                                                                 |
| Max Speed                 | 5 mph                                                                                                                                                                                    |
| Range Per Charge          | up to 20 mi.                                                                                                                                                                             |
| Motor                     | DC 24V / 300W x 2pcs                                                                                                                                                                     |
| Controller                | LiNX, PG VR2                                                                                                                                                                             |
| Battery Size (2 required) | 50Ah, 22NF AGM                                                                                                                                                                           |
| Charger                   | 5A off-board                                                                                                                                                                             |
| Gradient                  | 8°                                                                                                                                                                                       |
| Front Wheels              | 6" - solid                                                                                                                                                                               |
| Drive Wheels              | 12" - flat-free                                                                                                                                                                          |
| Rear Wheels               | 6" - solid                                                                                                                                                                               |
| HCPCS Codes               | K0824 / K0825<br>K0830 / K0831                                                                                                                                                           |
| Color Options             |  Red<br> Black |

Features MK batteries



merits



 **Enhancing life through mobility.**

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0425