


LATERAL PELVIC/ THIGH SUPPORT

PC's pelvic/thigh supports are universal and can accommodate all users with varying muscle tone and positioning needs.


Its adjustable hardware and multiple pad options can be configured for leg adduction and/or pelvic rotation.

This versatile solution is engineered for optimal lower extremity alignment, promoting pelvic stability for effective postural control.




Also comes as a removable hip guard


Key Features / Clinical Justification



Correct or prevent hip abduction/external rotation and secure pelvic stability



Varied offset range to accommodate windswept position



Support lower extremity alignment to promote better pressure distribution

HPCPS E1034 for lateral hardware, E0956 for positioning pads

Specification

- 1.25" extension length
- ~1.25" offset
- 300 lb. weight capacity
- Meets ISO16840-3 requirements

Available Pad Sizes

3"x5"
4"x6"
4"x8"
4"x12"