

ROADSTER

merits

RevVA9

SCOOTERS



ROADSTER 3

HCPCS Code

K0800

Weight Capacity: 300 lbs

Range: Up to 6 miles

Heaviest Piece: 29 lbs

FEATURES

- Extra charging port underneath the dashboard
- Detachable battery pack allows you to charge on the go
- Three wheels provide a tighter turning radius for better indoor maneuverability
- Designed to easily disassemble into five light components for transport and storage
- Adjustable tiller angle for ergonomic driving position and easy transferring onto and off of the unit
- 18" width deluxe stadium seat



ROADSTER 4

HCPCS Code

K0800

Weight Capacity: 300 lbs

Range: Up to 6 miles

Heaviest Piece: 36 lbs

FEATURES

- Extra charging port underneath the dashboard
- Detachable battery pack allows you to charge on the go
- Four wheels provide more stability
- Designed to easily disassemble into five light components for transport and storage
- Adjustable tiller angle for ergonomic driving position and easy transferring onto and off of the unit
- 18" width deluxe stadium seat

10/21/2021

merits

Merits Health Products Inc.

4245 Evans Avenue, Fort Myers FL, 33901

Phone: 1-800-963-7487 Fax: 239-772-3252

sales@meritsusa.com

	Roadster 3	Roadster 4
Model Number	S731A	S741A
Weight Capacity	300 lbs	300 lbs
Seat	18"x17"x15"	18"x17"x15"
Seat-to-Deck Height	16"-18"	16"-18"
Overall Length	40"	38"
Overall Width	21"	22"
Overall Height	36"	36"
Turning Radius	35.4"	43.3"
Ground Clearance	2"	2"
Weight of Heaviest Piece	29 lbs	36 lbs
Total Weight	99 lbs	106 lbs
Max Speed	4 mph	4mph
Per Charge Range	up to 6 miles	up to 6 miles
Batteries	12V/12Ah	12V/12Ah
Gradient	6°	6°
HCPCS Code	K0800	K0800
Warranty	Motors and electronics - 18 months, batteries - 1 year and frame - 5 years	

OPTIONS

- Oxygen Tank Holder
- Oxygen Tank Bag Holder
- Cane & Crutch Holder
- Walker Holder
- Foldable Rear Basket (Med & Large)
- Canopy
- Cover
- Mounting Bracket for Cup Holder
- Side Bag
- Safety Flag
- USB Charger
- Color Option: Red and Blue



Roadster 3



Merits Health Products



Meritshealthusa



Merits USA

merits

Merits Health Products Inc.

4245 Evans Avenue, Fort Myers FL, 33901

Phone: 1-800-963-7487 Fax: 239-772-3252

sales@meritsusa.com