

merits

# JUNIOR

K0822/K0823/K0815/50816

PORTABLE



# FEATURES

- Dynamic LiNX controller for easy programming and diagnostics
- 18" x 16" ergonomic, semi-reclining, high-back seat: up to 135°
- 9" knobby drive tires for superior traction
- Compact, light weight and easily disassembles into three pieces
- Interchangeable red and blue color shrouds

# SPECIFICATIONS

Weight Capacity	300 lbs.
Max Speed	4.5 mph
Per-Charge Range	12.5 miles
Batteries	12V/22Ah
Ground Clearance	2"
Turning Radius	23"
Overall Length	35"
Overall Width	23"
Drive Wheels	9"
Casters	7"
Seat Width	16" or 18"
Seat Depth	16"
Back Height	18"
Controller	Dynamic LinX
Total Weight	131 lbs
Weight of Heaviest Piece	61 lbs
<b>Warranty</b>	<b>18 months : Motor &amp; Electronics 5 Year/ Limited Lifetime Frame</b>

# OPTIONS

- Swing-away joystick
- USB charger
- Oxygen tank holder
- Foldable rear basket
- Cane & crutch holder
- Walker holder
- Weather cover
- Elevating leg rests

# MODELS

Junior - G2, Standard, Captain	P320-2S C
Junior - G2, Standard, Pan	P320-2S S
Junior - G1, Standard, Captain	P320-1S C
Junior - G1, Standard, Pan	P320-1S S



- Merits Health Products
- Meritshealthusa
- Merits USA

Merits Health Products Inc.  
 4245 Evans Avenue, Fort Myers FL, 33901  
 239-772-0579  
 sales@meritsusa.com