

PORTABLE POWER WHEELCHAIRS



EZ-GO

K0813/K0814

Weight Capacity: 250 lbs.
Range: Up to 9 miles

OPTIONS

Oxygen Tank Holder
USB charger
Cane and Crutch Holder
Walker Holder
Foldable Rear Basket
Cup Holder
Colors: red or blue

FEATURES

- Ergonomic mid-back seat for additional back support: 18W x 16D x 18H
- Quick connect feature for easy assembly and transportation
- Compact and lightweight design: overall weight 100 lbs.
- Turning radius of 23-inch for easy maneuverability in tight spaces
- Weight of the heaviest piece: 42 lbs.



EZ-GO DELUXE

RETAIL

Weight Capacity: 300 lbs.
Range: Up to 12 miles

OPTIONS

Oxygen tank holder
USB charger
Cane and crutch holder
Walker holder
Foldable rear basket
Cup holder
Colors: white, pink or turquoise

FEATURES

- Ergonomic mid-back seat for additional back support: 18"W x 16"D x 18"H
- Quick connect feature for easy assembly and transportation
- Compact and lightweight design with overall weight 113 lbs.
- 22 Ah batteries for extended range
- Turning radius of 23-inches for easy maneuverability in tight spaces
- Weight of the heaviest piece: 54 lbs.



Merits Health Products



Meritshealthusa



Merits USA

Merits Health Products Inc.

4245 Evans Avenue, Fort Myers FL, 33901

239-239-0579 | sales@meritsusa.com

merits



JUNIOR

K0815/K0816
K0822/K0823

Weight Capacity: 300 lbs.
Range: Up to 12.5 miles

OPTIONS

Swing-away joystick
USB charger
Oxygen tank holder
Foldable rear basket
Cane & crutch holder
Walker holder
Weather cover
Elevating legrests
Colors: red or blue

FEATURES

- Dynamic LiNX controller for easy programming and diagnostics
- 18 x 16-inch ergonomic, semi-reclining high-back seat: tilt up to 135°
- 9-inch knobby drive tires for superior traction
- 22 Ah batteries for extended range
- Compact, lightweight and easily disassembles into three pieces
- Weight of the heaviest piece: 61 lbs.



VISION CF

K0815/K0816
K0820/K0821/K0822/K0823

Weight Capacity: 300 lbs.
Range: Up to 18 miles

OPTIONS

Swing-away joystick
USB charger
Oxygen tank holder
Foldable rear basket
Cane & crutch holder
Walker holder
Weather cover
Elevating legrests

FEATURES

- 18 or 20-inch seat width with 18- 20-inch depth adjustability
- Ergonomic, semi-reclining high-back seat; up to 135°
- 9-inch drive tires, 6-inch casters
- Dynamic LiNX controller for convenient programming and diagnostics
- Premium AGM 34 Ah U1 batteries
- Weight of the heaviest piece:
53.5 lbs. (K0820/K0821)
63 lbs. (K0815/K0816/K0822/K0823)



TRAVEL-EASE

K0822

Weight Capacity: 300 lbs.
Range: Up to 18 miles

OPTIONS

Elevating legrests
USB charger
Oxygen tank holder
Seat cushion
Weather cover
Swing-away joystick

FEATURES

- 12.5-inch knobby drive tires for superior traction
- 4.8-inch ground clearance for easy obstacle clearance
- Height-adjustable flip-up armrests and swing-away footrests
- Foldable for improved storage and transportation
- Seat sizes available: 16, 18 or 20-inch
- Weight of the heaviest piece: 82 lbs.