

# VISION SUPER

HEAVY DUTY



merits

# FEATURES


- True mid-wheel drive: 20.7" turning radius
- Robust independent suspension can absorb up to 350 lbs on the footplate without tipping forward
- Premium, semi-reclining ergonomic seat, up to 135°
- 12" knobby drive tires for superior traction
- Dynamic LiNX & Shark controller for convenient Bluetooth programming and diagnostics
- 20" seat width with 18"-20" depth adjustability
- 22" or 24" seat widths with 20"-22" depth adjustability
- Red and blue color fenders included

# SPECIFICATIONS

Weight Capacity	450 lbs, 400 lbs (with power seat lift)
Maximum Speed	5 mph
Per-Charge Range	20 miles
Batteries	22NF (55Ah)
Ground Clearance	2.2"
Turning Radius	20.7"
Overall Length	35"
Overall Width	27"
Drive Wheels	12" flat-free
Casters	6" front, 6" rear
Suspension	TruStep
Seat-to-Floor Height	21"-23" (captain), 17.5"-19.5" (pan) Power seat lift adds 3"
Seat Width	20", 22" or 24"
Seat Depth	18", 19", 20", 21" or 22"
Backrest Height	20" or 22"
Power Seat Lift	10" scissor (P3274)
Controller	LiNX, PG VR2 (for power seat lift)
Warranty	18 Months : Motor and Electronics 5 Year/ Limited Lifetime Frame 1 Year Battery

# OPTIONS

- Unoccupied transit tie down kit
- 24"W x 20"-22"D x 22"H seat
- 10" scissor power seat lift
- Swing-away joystick
- USB Charger
- Elevating legrests
- Oxygen tank holder
- Foldable rear basket
- Walker holder
- Cane & crutch holder
- Color Options: Red OR Blue



# MODELS

Vision Super - Captain	P327-2HD C
Vision Super- Pan	P327-2HD S
Vision Super w/PSL- Captain	P3274-2SSE C
Vision Super w/PSL- Pan	P3274-2SSE S



- Merits Health Products
- Meritshealthusa
- Merits USA

Merits Health Products Inc.  
 4245 Evans Avenue, Fort Myers FL, 33901  
 239-772-0579 sales@meritsusa.com