

VISION

K0837

ULTRA HD



merits

© Merits Health Products, Inc. Rv3.22-01

FEATURES

- Mid-wheel drive design provides outstanding maneuverability
- Adjustable seat width 20"-24" and seat depth 16"-22" for custom fit
- Optimized seating stability with 4-post seat design
- Tool-free battery box provides easy maintenance
- Front and rear suspension offers all day comfort and easy obstacle clearance
- Fixed back angle: 90°, 95°, 100°, 105° and 110°
- Dual in-line motors for enhanced efficiency, torque, range, and performance
- Unoccupied transit tie down brackets included

SPECIFICATIONS

Weight Capacity	450 lbs
Maximum Speed	5 mph
Per-Charge Range	8 miles/15 miles
Batteries	34Ah/U1 or 55Ah/22NF
Controller	PG VR2-A2 60A
Ground Clearance	2.875"
Turning Radius	21"
Overall Length	35"/ 45" incl. Footplate
Overall Width	27"
Drive Wheels	12" x 3"
Casters	6" x 2"
Suspension	Front and Rear
Seat Width	20"-24", adjusted in 1" increments
Seat Depth	16"-22", adjusted in 1" increments
Backrest Height	21"-23"
Seat-to-Floor Height	17.5" or 18.5"
Overall Weight	288 lbs (without front riggings) 310 lbs (with Power AFP)
HCPSCS	K0837
Warranty	18 Months : Motor and Electronics 5 year/ limited lifetime frame 1 Year Battery

OPTIONS

- Power Articulating Foot Platform
- Elevating Legrests
- Power Elevating Legrests
- Trunk Supports: Available in 3 Sizes
- Abductor: Available in 2 Sizes
- Swing-away Joystick Mounting Bracket
- Battery Upgrade - 22NF
- Color Options:



Merits Health Products Inc.

4245 Evans Avenue, Fort Myers FL, 33901

Phone: 1-800-963-7487 Fax: 239-772-3252

Email: sales@meritsusa.com

merits