

VISION

ULTRA

K0822/K0823/K0835



merits

Merits Health Products, Inc. Rv6.22-01

FEATURES

- Large 14" knobby drive tires
- Tension-adjustable front and rear suspension
- Unoccupied tie down restraints
- Multiple seating options, power tilt, rehab seat frame, captain or solid seat pan
- Multiple seat-to-floor height options: 16.5"-19.5"
- Width and depth-adjustable seat; 16"-22" in 1" increments
- Red, white, blue and charcoal gray color options

SPECIFICATIONS

Weight Capacity	300 lbs
Maximum Speed	5 mph
Per-Charge Range	12 miles (19 miles w/22NF)
Batteries	34Ah/U1 or 55Ah/22NF
Ground Clearance	2.4"
Turning Radius	26"
Overall Length	36"
Overall Width	24"
Drive Wheels	12" x 3"
Casters	6" x 2"
Suspension	Front and Rear
Seat-to-Floor Height	Power tilt or pan: 16.5"-19.5", Captain: 20"-23"
Seat Width	16"-22", adjusted in 1" increments
Seat Depth	16"-22", adjusted in 1" increments
Backrest Height	21"-23"
Overall Weight	315 lbs (22"W Captain Seat, 22NF Batteries) 332 lbs (Power Tilt, Power AFP, 22NF Batteries)
Controller	PG VR2-A2 60A
Warranty	18 Months: Motor and Electronics
	5 year/ limited lifetime frame
	1 Year Battery


OPTIONS

- 22NF battery upgrade
- Power articulating foot platform
- Elevating leg rests or swing-away footrests
- Swing-away joystick
- USB charger
- Weather cover



MODELS

Vision Ultra - G2, Standard, Captain	P325-2S C
Vision Ultra - G2, Standard, Pan	P325- 2S S
Vision Ultra - G2, Standard, Single Power	P325-2SSP S

 Merits Health Products

 Meritshealthusa

 Merits USA

Merits Health Products Inc.

4245 Evans Avenue, Fort Myers FL, 33901

239-239-0579 | sales@meritsusa.com

merits