

Roadster⁴



S741RS4-SBB / Galaxy Blue
S741RS4-ARB / Cherry Red

SPECIFICATIONS

Product #	S741RS4
Weight Capacity	300 lbs.
Range	up to 12 mi.
Seat Width	18"
Seat Depth	18"
Backrest Height	14.5" from seat
Seat-to-Ground Height	23.5"
Seat -to-Deck Height	18.5"
Overall Length	40.5"
Overall Width	22"
Overall Height	36.6"
Turning Radius	47"
Ground Clearance	2.5"
Total Weight	107 lbs.
Weight of Heaviest Piece	36 lbs.
Max Speed	4 mph
Basket Size	12.5Wx10Hx6D
Motor	24V 220W
Battery Size (2 required)	12AH
Charger	29V2AH
Gradient	6°
Front Wheel (s)	7"
Rear Wheels	8"
HCPCS Code	na
Standard Colors	Cherry Red and Galaxy Blue



- Maximum user weight 300 lbs.
- Maximum speed 4 mph.
- Drive range up to 12 miles.
- 5-piece design for easy assembly and disassembly. (no tools required)
- Weight of heaviest piece - 36 lbs.
- Headlight for visibility.
- Height adjustable seat with swivel.
- Convenient storage pocket.

MSRP **\$2,110**



4245 Evans Ave., Fort Myers, FL 33901

239-772-0579

W www.meritsusa.com E sales@meritsusa.com