

# YOGA

# SCOOTER



Foldable scooter includes convenient carrying case.



**merits**

Merits Health Products Inc.  
4245 Evans Avenue | Fort Myers FL | 33901  
239-772-0579

# FEATURES

- Strong light-weight aluminum alloy frame for easy folding and carrying
- Small suitcase with wheels for easy transport and storage
- Easy to connect lithium-ion battery
- Height and angle adjustable steering tiller to fit individual user demand
- Patented twin-front-wheels with spring suspension system offers excellent stability

# ORDER FORM

Provider Name \_\_\_\_\_

Client ID \_\_\_\_\_

Acct # \_\_\_\_\_ Date \_\_\_\_\_

Ship to \_\_\_\_\_

City \_\_\_\_\_

State \_\_\_\_\_ Zip \_\_\_\_\_

Phone \_\_\_\_\_

Email \_\_\_\_\_

PRODUCT	PART#	MSRP
Yoga with Case	S542	\$3964.00

Quantity \_\_\_\_\_

# SPECIFICATIONS

Weight Capacity	250 lbs
Max Speed	4 mph
Per-Charge Range	up to 11.25 miles
Turning Radius	37.8"
Ground Clearance	3.7"
Motor	DC 2v,80w 340w (MAX)
Controller	PG S - Drive 45A
Battery	24V / 11.5Ah Lithium-ion
Charger	4A off-board
Gradient	6°
Front Wheels	6" (153x40) foam filled tires
Rear Wheels	7"(178x45) foam filled tires
Brakes	Intelligent, Regenerative, Electromagnetic
Warranty	18 Months : Motor & Electronics 1 Year Battery



Merits Health Products Inc.  
4245 Evans Avenue | Fort Myers FL | 33901  
239-772-0579

**merits**