

YOU ASKED FOR IT.

WE GOT IT.

ULTRA 3

TILT | ELEVATE
ARTICULATING FOOT PLATFORM

K0835



Group 2
Single Power
3 Actuators

merits

ULTRA 3

TILT | ELEVATE | AFP
K0835

- Power tilt 0° to 50°
- Power seat lift with 12" telescoping tower
- Articulating foot platform with simultaneous extension and articulation
- Tension adjustable front and rear suspension for comfort and easy obstacle clearance
- Multiple seat-to-floor height options 16.5"-19.5"
- Adjustable seat width and depth 16" - 22"



Dimensions	
Seat Width	16"-22" adj. in 1" increments
Seat Depth	16"-22" adj. in 1" increments
Backrest Height	21"-23" (18"-28" ProForm)
Seat-to-Floor Height (no cushion)	16.5"-19.5"
Leg Length (AFP)	10"-18"
Overall Length	36"
Overall Width	24"
Overall Height	40"
Lift	12" telescoping
Technical Specs	
Weight Capacity	300 lbs
Maximum Speed	5 mph
Per-Charge Range	up to 12 miles (19 miles w/22NF)
Turning Radius	26" including footplates
Ground Clearance	2.5"
Motor	DC 24V, 240W
Controller	PG 24W VR2-A2 60A
Batteries	34Ah/U1 or 55Ah/22NF
Charger	6A off-board
Gradient	6°
Caster	6" x 2"
Drive Wheel	14"x3"
Base	150 lbs
Overall Weight	384 lbs (Elevate, Tilt, AFP, 22NF Batts)
Warranty	18 Months: Motor & Electronics
	5 year/ limited lifetime frame
	1 Year Battery

ULTRA

TILT | CAPTAIN SEAT
K0836



ULTRA X

330 LB CAPACITY
K0837



239-772-0579

NOTE: Elevate is not available on Ultra Captain/Tilt or UltraX.

MERITS HEALTH PRODUCTS | 4245 EVANS AVE | FORT MYERS FL 33901 | WWW.MERITSUSA.COM