

merits

DUALER

K0830/K0831

FULL SIZE



FEATURES


- Unit changes from rear to front-wheel drive with the swivel of a seat
- 5-inch power seat lift
- 16 or 18-inch ergonomic semi-reclining seats; up to 135°
- 10-inch knobby drive tires for superior traction
- Dual footplates for comfort and convenience
- Available in Red

SPECIFICATIONS

Weight Capacity	300 lbs
Max Speed	5 mph
Per-Charge Range	up to 18 miles
Batteries	U1-34Ah x 2
Ground Clearance	2"
Turning Radius	21"
Overall Length	38"
Overall Width	24"
Drive Wheels	10" foam filled tires
Rear Casters	8" foam filled tires
Seat-to-Floor Height	21.5-26"
Seat Width	16", 18"
Seat Depth	16" - 18"
Back Height	18" or 20"
Controller	PG VR2
Power Seat Lift	5"- single post
Warranty	18 Months : Motor & Electronics 5 Year/ Limited Lifetime Frame 1 Year/ Batteries

MODELS

Dualer - FWD/RWD, Captain	P312-2C
Dualer - FWD/RWD, Pan	P312-2S

 Merits Health Products

 Meritshealthusa

 Merits USA

OPTIONS

- Swing-away joystick
- USB charger
- Oxygen tank holder
- Foldable rear basket
- Cane & crutch holder
- Walker holder
- Weather cover
- Elevating legrests



Merits Health Products Inc.

4245 Evans Avenue, Fort Myers FL, 33901

239-239-0579 | sales@meritsusa.com

merits