

ULTRA CST

CAPTAIN SEAT | TILT

K0836



merits

ULTRA CST

CAPTAIN SEAT | TILT

K0836

- Manual recline adjustment from 90° to 105°
- **Power tilt** from 0° to 50°
- The center of gravity changes with the seat angle (keeps the center of gravity above the drive wheels)
- Adjustable height removable armrests
- Seat width 18-inch or 20-inch
- Adjustable seat depth (18-, 19-, or 20-inch for customized fit)
- Full front and rear suspension

OPTIONAL

- Power footplate (AFP)



Dimensions	
Seat Width	18" or 20"
Seat Depth	18"-20" adj. in 1" increments
Backrest Height	20"
Seat-to-Floor Height (no cushion)	21.5-22.5
Overall Length	36"
Overall Width	24"
Overall Height	40"
Technical Specs	
Weight Capacity	300 lbs
Maximum Speed	up to 5 mph
Per-Charge Range	up to 11 miles (19 miles w/22NF)
Turning Radius	28" including footplates
Ground Clearance	2.5"
Motor	DC 24V, 240W
Controller	PG 24W VR2-A2 60A
Batteries	24v/22NF*
Charger	6A off-board
Gradient	6°
Caster	6" x 2"
Drive Wheel	14"x3"
Base	128 lbs
Seat Weight	99 lbs
Overall Weight	310 lbs (with Batteries)
Warranty	18 Months: Motor & Electronics
	5 year/ limited lifetime frame
	1 Year Battery

*Recommended for best performance

AVAILABLE IN RED



merits *AviD* *Pilot*

239-772-0579

sales@meritsusa.com