

# AXCEL NARROW



**22.5**  
**INCHES WIDE**

*The Industry's First 22.5-inch Wide  
Front Wheel Drive Wheelchair*

**Avid** Rehab

[www.meritususa.com](http://www.meritususa.com)

# AXCEL NARROW

With its 14" front drive wheels and its narrow pivot points, the AXCEL Narrow gives you the outdoor performance of FWD and the maneuverability of a MWD!

- 22.5-inch wide front wheel drive chair is easy to maneuver tight spots like narrow hallways and door frames.
- Tight 90° turns
- Navigating over uneven terrain is a snap with FWD because your 14" drive tires maintain contact.
- High-performance, independent hydraulic suspension
- Large 9-inch rear casters to navigate obstacles and rough terrain
- Bright LED front and rear lights
- Tool-free rear battery access
- Easy access for electronics
- Power recline adjustment 85°-160°
- 12" power seat lift option



Seat Width	14" - 20"
Seat Depth (1.5" PC backrest cushion)	14" - 22"
Seat-to-Floor Height (w/o cushion)	17.5" - 19.5"
Leg Length (AFP)	10" - 18"

Overall Length	39.6"
Overall Width	22.5"
Overall Height	41.5"

Weight Capacity	300 lbs.
Max. Speed	Up to 5.5mph
Range	12 mi.
Turning Radius (w/o footplates)	26"
Ground Clearance	3"
Motor	5300rpm / 222W(1:29)
Controller	R-Net 90A
Battery	22NF 50AHx 2 pcs
Charger	6A off-board
Gradient	7.5°
Caster	9"
Drive Wheel	14" x 3"
Base	143 lbs.
Battery/ea	34 lbs.
Seat (Tilt & Recline)	122.5lbs.
3 Actuators	385 lbs

*weight may vary based on selected options*



RUBY | SAPPHIRE | CAMO | DOVE | NOIR

**K0861/K0856**



**239-772-0579**

© Merits Health Products, Inc 2024