

EZ-GO



 **merits**

Enhancing life through mobility.

EZ-GO

- ✦ Ergonomic mid-back seat for additional back support (18" x 17" x 18")
- ✦ Higher weight capacity of 300 lbs.
- ✦ Quick connect feature for easy assembly and transportation
- ✦ Compact and light weight design
- ✦ 22 AH batteries for extended range
- ✦ Easy maneuverability in tight spaces with a turning radius of 23"

Product #	P321B
Weight Capacity	300 lbs
Drivetrain	RWD
Seat Type	Stadium
Seat Width	18"
Seat Depth	18"
Backrest Height	17"
Seat-to-Floor Height	18" - 20"
Elevating Seat	n/a
Overall Length	27"
Overall Width	23"
Overall Height	43"
Turning Radius	23"
Ground Clearance	1.8"
Total Weight	113 lbs
Wt of Heaviest Pc	54 lbs
Max Speed	4.5 mph
Per Charge Range	up to 12 miles
Motor	DC 24V, 90W
Controller	P&G NR2
Battery Size (2)	12V/22AH
Charger	1.5A off-board
Gradient	6°
Front Wheels	7" solid
Drive Wheels	8" foam filled
Rear Wheels	3" anti-tippers
Group	Group: 1
HCPCS Code	K0814



Black



Red

 **merits**

merits *AviD* *Pilot*

239-772-0579