

Vision Ultra HD Order Form


PROVIDER INFORMATION

Provider Name _____
 ATP _____
 Client ID _____
 Account Number _____
 Phone _____
 Fax _____
 Email _____
 Date _____
 PO# _____
 Ship To Address _____
 City _____
 State _____
 Zip _____



CLIENT INFORMATION

A. Hip Width _____
 B. Seat Depth _____
 C. Lower Leg Length _____
 D. Seat to Top of Head _____
 E. Seat to Top of Shoulder _____
 F. Seat to Axilla (L,R) Armpit _____
 G. Seat to Arm _____
 H. Chest Width _____
 I. Shoulder Width _____
 J. Client Height _____
 K. Client Weight _____



DME Providers are responsible for determining appropriate billing codes when submitting claims for reimbursement. HCPCS codes should not be used as legal or billing advice.

Base	Part Number	HCPCS	MSRP
<i>Group 2 Power Base, with 2 pole 200W motors, 450lbs weight capacity</i>			
<input type="checkbox"/> Vision Ultra HD Power Base	30200671	K0837	\$8,249.00
<i>^^^Unoccupied transit tie downs are included^^^</i>			

Color	Part Number	HCPCS	MSRP
<input type="checkbox"/> Candy Apple Red	30500858	-	No Charge
<input type="checkbox"/> White	30500859	-	No Charge
<input type="checkbox"/> Blue	30500860	-	No Charge
<input type="checkbox"/> Charcoal Gray	30500861	-	No Charge

Batteries	Part Number	HCPCS	MSRP
<input type="checkbox"/> 12V, 34 Ah, U1, AGM, BWG (x2) Each	SCPU1	E2365	\$146.30
<input type="checkbox"/> 12V, 55 Ah, 22NF, AGM, BWG (x2) Each	BW12550-IT	E2361	\$280.50
<input type="checkbox"/> 12V, 55 Ah, 22NF, Gel, MK (x2) Each	M22NF	E2361	\$324.50
<i>^^^ We recommend using 22NF batteries^^^</i>			

Electronics	Part Number	HCPCS	MSRP
<input type="checkbox"/> Through Toggle Box , Adjustable	33700589	-	No Charge
<input type="checkbox"/> Through VR2 - One Actuator	P3230001-NI	E2310	\$1,100.00
<input type="checkbox"/> Through VR2 - Two + Actuators	33500524-NI	E2311	\$2,895.00

Joystick Mounting Bracket	Part Number	HCPCS	MSRP
<input type="checkbox"/> In-line Joystick Mounting Bracket <input type="checkbox"/> Right <input type="checkbox"/> Left	33500883	-	\$275.00
<input type="checkbox"/> Swing-away Mounting Bracket (Right)	39922019-NI	E1028	\$275.00
<input type="checkbox"/> Swing-away Mounting Bracket (Left)	39922020-NI	E1028	\$275.00
<input type="checkbox"/> Retract4 Swing-away Joystick Mount <input type="checkbox"/> Right <input type="checkbox"/> Left	70-JS460/522	E1028	\$350.00

Joystick Knob Options	Part Number	HCPCS	MSRP
<input type="checkbox"/> Joystick Knob, 80MM Stick	SA72435	E2323	\$80.00
<input type="checkbox"/> Joystick Knob, Mushroom	SA72439	E2323	\$88.00
<input type="checkbox"/> Joystick Knob, 38MM Ball	SA72441	E2323	\$64.00
<input type="checkbox"/> Joystick Knob, Chin	SA75165	E2323	\$64.00
<input type="checkbox"/> Joystick Knob, T-Bar	SA80941	E2323	\$132.00
<input type="checkbox"/> JoystickKnob, Soft Ball	SA81095	E2323	\$100.00

Power Tilt Seating	Part Number	HCPCS	MSRP
<input type="checkbox"/> Power Tilt Seat/Angle-Adjustable Back	32200657	E1002	\$6,049.00

Seat to Floor Height Options					
		<input type="checkbox"/> 17.5"			<input type="checkbox"/> 18.5"
Dimensions for Seating					
Width			Depth		
Width measured in 1" Increments	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	Depth measured in 1" Increments	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"
	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"		<input type="checkbox"/> 18"	<input type="checkbox"/> 19"
	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"		<input type="checkbox"/> 20"	<input type="checkbox"/> 21"
	<input type="checkbox"/> 24"			<input type="checkbox"/> 22"	

Privacy Flap Options	Part Number	HCPCS	MSRP
<input type="checkbox"/> Privacy Flap for Seat Widths of 20" - 24"	43780280	-	\$28.00

PC Gel		Part Number	HCPCS Code	MSRP
<input type="checkbox"/>	PC Gel 18"x16"	G1816	E2603	\$299.00
<input type="checkbox"/>	PC Gel 18"x18"	G1818	E2603	\$299.00
<input type="checkbox"/>	PC Gel 18"x20"	G1820	E2603	\$299.00
<input type="checkbox"/>	PC Gel 20"x16"	G2016	E2603	\$299.00
<input type="checkbox"/>	PC Gel 20"x18"	G2018	E2603	\$299.00
<input type="checkbox"/>	PC Gel 20"x20"	G2020	E2603	\$299.00
<input type="checkbox"/>	PC Gel 22"x18"	G2218	E2604	\$359.00
<input type="checkbox"/>	PC Gel 22"x20"	G2220	E2604	\$359.00
<input type="checkbox"/>	PC Gel 22"x22"	G2222	E2604	\$359.00
<input type="checkbox"/>	PC Gel 22"x24"	G2224	E2604	\$359.00
<input type="checkbox"/>	PC Gel 24"x18"	G2418	E2604	\$359.00
<input type="checkbox"/>	PC Gel 24"x20"	G2420	E2604	\$359.00

PC Lite		Part Number	HCPCS Code	MSRP
<input type="checkbox"/>	PC Lite 18"x16"	L1816	E2601	\$89.00
<input type="checkbox"/>	PC Lite 18"x18"	L1818	E2601	\$89.00
<input type="checkbox"/>	PC Lite 18"x20"	L1820	E2601	\$89.00
<input type="checkbox"/>	PC Lite 20"x16"	L2016	E2601	\$89.00
<input type="checkbox"/>	PC Lite 20"x18"	L2018	E2601	\$89.00
<input type="checkbox"/>	PC Lite 20"x20"	L2020	E2601	\$89.00
<input type="checkbox"/>	PC Lite 22"x18"	L2218	E2602	\$169.00
<input type="checkbox"/>	PC Lite 22"x20"	L2220	E2602	\$169.00
<input type="checkbox"/>	PC Lite 22"x22"	L2222	E2602	\$169.00
<input type="checkbox"/>	PC Lite 22"x24"	L2224	E2602	\$169.00
<input type="checkbox"/>	PC Lite 24"x18"	L2418	E2602	\$169.00
<input type="checkbox"/>	PC Lite 24"x20"	L2420	E2602	\$169.00

Seat Belt		Part Number	HCPCS	MSRP
<input type="checkbox"/>	67" Seat Belt (standard)	43700319-NI	-	No Charge
<input type="checkbox"/>	84" Seat Belt/ Axle Mount	SB225	E0978	\$80.00

PC Deep		Part Number	HCPCS Code	MSRP
<input type="checkbox"/>	PC Deep 18"x16"	D1816	E2607	\$389.00
<input type="checkbox"/>	PC Deep 18"x18"	D1818	E2607	\$389.00
<input type="checkbox"/>	PC Deep 18"x20"	D1820	E2607	\$389.00
<input type="checkbox"/>	PC Deep 20"x16"	D2016	E2607	\$389.00
<input type="checkbox"/>	PC Deep 20"x18"	D2018	E2607	\$389.00
<input type="checkbox"/>	PC Deep 20"x20"	D2020	E2607	\$389.00
<input type="checkbox"/>	PC Deep 22"x18"	D2218	E2608	\$469.00
<input type="checkbox"/>	PC Deep 22"x20"	D2220	E2608	\$469.00
<input type="checkbox"/>	PC Deep 22"x22"	D2222	E2608	\$469.00
<input type="checkbox"/>	PC Deep 22"x24"	D2224	E2608	\$469.00
<input type="checkbox"/>	PC Deep 24"x18"	D2418	E2608	\$469.00
<input type="checkbox"/>	PC Deep 24"x20"	D2420	E2608	\$469.00

Accessories		Part Number	HCPCS	MSRP
<input type="checkbox"/>	USB Charger	61230167	-	\$45.00
<input type="checkbox"/>	Oxygen Tank Holder - Bag	34100241	-	\$179.00

Backrests		Part Number	HCPCS	MSRP
<input type="checkbox"/>	Omit back Pan	-	-	No Charge
<input type="checkbox"/>	18"W x 17"H, Adjustable Height Range 21.5"-23"	43780112/43780113	-	No Charge
<input type="checkbox"/>	20"W x 17"H Back Pan, Adjustable Height Range 21.5"-23"	43780128/43780130	-	\$164.00
<input type="checkbox"/>	22"W x 17"H Back Pan, Adjustable Height Range 21.5"-23"	43780193/43780194	-	\$164.00
<input type="checkbox"/>	24"W x 17"H Back Pan, Adjustable Height Range 21.5"-23"	43780192/43780195	-	\$164.00
<input type="checkbox"/>	Adapter for comfort company backrest	39900215	-	No Charge

*** Not for use with Avid backrest or headrest***

<input type="checkbox"/>	Add Endomorphic Back Kit -	39910068		
--------------------------	----------------------------	----------	--	--

Angle-adjustable back

90° 95° 100° 105° 110°

Precision Comfort ProForm Backrest		Part Number	HCPCS	MSRP
<input type="checkbox"/>	PC PF Planar Backrest 18"W x 18"H	PF1818	E2620	\$675.00
<input type="checkbox"/>	PC PF Planar Backrest 20"W x 18"H	PF2018	E2620	\$675.00
<input type="checkbox"/>	PC PF Planar Backrest 22"W x 18"H	PF2218	E2621	\$755.00
<input type="checkbox"/>	PC PF Planar Backrest 24"W x 18"H	PF2418	E2621	\$755.00
<input type="checkbox"/>	PC PF Planar Backrest 18"W x 22"H	PF1822	E2621	\$675.00
<input type="checkbox"/>	PC PF Planar Backrest 20"W x 22"H	PF2022	E2620	\$675.00
<input type="checkbox"/>	PC PF Planar Backrest 22"W x 22"H	PF2222	E2621	\$755.00
<input type="checkbox"/>	PC PF Planar Backrest 24"W x 22"H	PF2422	E2621	\$755.00

Upper Extremity Positioning Components		Part Number	HCPCS	MSRP
<input type="checkbox"/>	Swing-away Trunk Support Bracket Pair	39900179	E1028	\$549.00
<input type="checkbox"/>	Trunk Support Pad- Medium 5" x 4" Pair	39900177	E0956	\$87.00
<input type="checkbox"/>	Trunk Support Pad- Small 4" x 3" Pair	39900178	E0956	\$87.00

Headrest - REQUIRED		Part Number	HCPCS	MSRP
<input type="checkbox"/>	Omit Headrest	-	-	No Charge
<input type="checkbox"/>	PC ProForm Midline Headrest w/ 10" Pad Kit	32100388/32000640	E1028/E0955	\$410.00
<input type="checkbox"/>	PC ProForm Multi-Axis Headrest w/ 10" Pad Kit	32100372/32000640	E1028/E0955	\$430.00

Armrests	Part Number	HCPCS	MSRP
<input type="checkbox"/> Single Post, Height-Adjustable, Flip-Back Armrests	39910060/61	-	No Charge
Armrest Options	Part Number	HCPCS	MSRP
<input type="checkbox"/> Side Guard	34100277	-	\$50.00
Gel Ovations Arm Accessories	Part Number	HCPCS	MSRP
<input type="checkbox"/> Gel Ovations Arm Pads 3.5"x12" Pair	312G	E2209	\$140.00
<input type="checkbox"/> Gel Ovations Arm Pads 3.5"x14" Pair	314G	E2209	\$170.00
<input type="checkbox"/> Gel Ovations Arm Pad 2" X 12" Pair	212G	E2209	\$115.00
<input type="checkbox"/> Gel Ovations Arm Pad 2" X 14" Pair	214G	E2209	\$130.00
<input type="checkbox"/> Gel Ovations Universal Elbow Stop Bracket and Flat 3"x4" Pad Each	<input type="checkbox"/> Right <input type="checkbox"/> Left UESF1	E2209	\$107.00
Gel Ovations Arm Troughs	Part Number	HCPCS	MSRP
<i>Arm trough must be ordered with 3.5" and 2" Gel Armrests</i>			
<input type="checkbox"/> Gel Ovations 12" First Trough	<input type="checkbox"/> Right <input type="checkbox"/> Left FTS-312G	E2209	\$360.00
<input type="checkbox"/> Gel Ovations 12" Second Trough	<input type="checkbox"/> Right <input type="checkbox"/> Left STS-312G	E2209	\$360.00
<input type="checkbox"/> Gel Ovations 12" First & Second Trough	<input type="checkbox"/> Right <input type="checkbox"/> Left BTS-312G	E2209	\$580.00
<input type="checkbox"/> Gel Ovations 12" First & Second Trough w Elbow St	<input type="checkbox"/> Right <input type="checkbox"/> Left BTS-312G-UESF1	E2209	\$687.00
<input type="checkbox"/> Gel Ovations 14" First Trough	<input type="checkbox"/> Right <input type="checkbox"/> Left FTS-314G	E2209	\$415.00
<input type="checkbox"/> Gel Ovations 14" Second Trough	<input type="checkbox"/> Right <input type="checkbox"/> Left STS-314G	E2209	\$415.00
<input type="checkbox"/> Gel Ovations 14" First & Second Trough	<input type="checkbox"/> Right <input type="checkbox"/> Left BTS-314G	E2209	\$660.00
<input type="checkbox"/> Gel Ovations 14" First & Second Trough w Elbow St	<input type="checkbox"/> Right <input type="checkbox"/> Left BTS-314G-UESF1	E2209	\$767.00

Footrest Options	Part Number	HCPCS	MSRP
Foot Platform and Legrest Options	Part Number	HCPCS	MSRP
<input type="checkbox"/> Omit Front Riggings	-	-	No Charge
<input type="checkbox"/> Power Articulating Foot Platform - 14"-16" Lower Leg Length	33300472	E1012	\$2,398.00
> <input type="checkbox"/> Large Footplate - 13.5" x 11.5"	-	-	No Charge
> <input type="checkbox"/> Short Bracket 13"-15" Lower Leg Length	39933059	-	No Charge
<i>****Only compatible with the large footplate. It is not compatible with the split footplates / Only applicable with 18.5" seat-to-floor height****</i>			
> <input type="checkbox"/> Split Footplate	-	-	-
> <input type="checkbox"/> 5" x 6" Pair	33300481	K0040	\$145.00
> <input type="checkbox"/> 3" Extension Pair	33300489	E2209	\$50.00
> <input type="checkbox"/> Side Stop Plate for 5" x 6" Pair	33300495	E2209	\$25.00
> <input type="checkbox"/> 5" x 11" Pair	33300473	K0040	\$160.00
> <input type="checkbox"/> 3" Extension Pair	33300489	E2209	\$50.00
> <input type="checkbox"/> Side Stop Plate for 5" x 11" Pair	33300493	E2209	\$25.00
<input type="checkbox"/> Power Elevating Legrests Pair; 5"W x 6"D footplate	33400296	E1010	\$2,398.00
Manual Footrest Options	Part Number	HCPCS	MSRP
<input type="checkbox"/> Seat-mount Foot Platform w/Calf Pad -13"-15" Lower Leg Length, 13.5" x 11.5" Footplate	33300398	-	No Charge
<input type="checkbox"/> Swing-away Footrest - 13" - 17" ext. range (adj. in 0.5" increments) 5.5" x 6" Footplates w/ Heel Loops, 70° Angle, Pair	33300424	-	No Charge
<input type="checkbox"/> Manual Elevating Legrests- 13"-18" Adjustment Range; 5"W x 6"D footplate	33400297	E0990	\$499.00
<input type="checkbox"/> Sta-Rite Residual Limb Support 6" x 7.5" Gel Pad for Flat Mount	<input type="checkbox"/> Right <input type="checkbox"/> Left SRSSF	E1020	\$350.00
<input type="checkbox"/> Sta-Rite Residual Limb Support Swing-away Mounting Hardware	<input type="checkbox"/> Right <input type="checkbox"/> Left SRSSFHO	E1028	\$330.00
Lower Extremity Positioning Components	Part Number	HCPCS	MSRP
<input type="checkbox"/> Gel Ovations Foot Platform Gel Cover 11" x 9"	GFPCS	E2209	\$150.00
<input type="checkbox"/> Gel Ovations Gel and Foam Calf Pads - 8" x 6" x 1.5" Pair	GCP	E2209	\$245.00
<input type="checkbox"/> Sta-Rite Axis Adjustable Hip Guides 3.5" x 12" Gel Pad	<input type="checkbox"/> Right <input type="checkbox"/> Left SRMHG12	E0956	\$330.00
<input type="checkbox"/> Hip Guide Swing-Away Mounting Hardware	SRMHG12-NI	E1028	\$248.00
<input type="checkbox"/> Sta-Rite Adductor Gel Pad, 3" x 4"	<input type="checkbox"/> Right <input type="checkbox"/> Left SRASF1	E0957	\$385.00
<input type="checkbox"/> Adductor Bracket Swing-Away Mounting Hardware	SRASF1-NI	E1028	\$248.00
<input type="checkbox"/> Abductor Pad - Medium - 6" x 4" x 3"	39900174	E0957	\$59.00
<input type="checkbox"/> Abductor Pad - Small - 5" x 3" x 3"	39900175	E0957	\$57.00
<input type="checkbox"/> Flip-Down and Removable Abductor Bracket	39900176	E1028	\$164.00