



S941 4-Wheel Full Suspension Electric Scooter

S941

4-Wheel Full Suspension Electric Scooter



Features

- Brand new sporty streamline design.
- Full suspension system with shock absorbers delivers a smooth and comfortable ride.
- Modern LED lighting system delivers a brilliant light with less power consumption.
- Digital LCD dashboard display: time display, battery gauge, temperature, maintenance messages, speed display and odometer.
- Extra storage space beneath tiller.
- Automatic safety turning half speed control by sensor.



Modern LED lighting system delivers brilliant light with less power consumption.

Dimension(cm / in)

seat width	seat depth	seat height from ground
46/18	40/16	63.5/25

seat height from deck	back height	wheel base
43-48/17-19	46/18	97/38

length	width	height
130/51	69/27	125/49



More space for storage



Digital LCD dashboard display

Technical Specs.

max. speed up to	12kph / 7.5mph ; 15kph / 9.5mph
range up to	40km / 25mi 167cm / 65.7" 9cm /
turn radius	3.5"
ground clearance	DC24V,400W(Rating)/1750W(Max)
motor	PG S-Drive 120 ; Dynamic Rhino2-
controller	120A 12V 50AH x 2 5A off-board up
battery	to 10 ° 12" 12" intelligent,
charger	regenerative,
gradient	electromagnetic brakes.
front wheel	
rear wheel	
brake	

Weight(kg / lb)

base	battery/ea	seat	tll w/batt.	capacity
69/152	17/37	18/40	118/260	159/350



730 N.E. 19th Place, Cape Coral, FL 33909
 TEL:1-800-963-7487, 1-239-772-0579
 SALES FAX:1-239-574-2661 MAIN FAX:1-239-772-3252
 www.meritsusa.com