

ULTRA X

K0837



merits

Merits Health Products, Inc. 103023


FEATURES

- 330 lb. capacity
- Large 14" knobby drive tires
- Tension-adjustable front and rear suspension
- Unoccupied tie down restraints
- Multiple seat-to-floor height options: 16.5"-19.5"
- Width and depth-adjustable seat; 16"-22" in 1" increments
- Red, white, blue and charcoal gray color options


SPECIFICATIONS

Weight Capacity	330 lbs
Maximum Speed	5 mph
Per-Charge Range	up to 26 miles
Batteries	34Ah/U1 or 55Ah/22NF
Ground Clearance	2.5"
Turning Radius	26" (including footplates)
Overall Length	36"
Overall Width	24"
Drive Wheels	14" x 3"
Casters	6" x 2"
Suspension	Front and Rear
Seat-to-Floor Height	Power tilt 16"-19"
Seat Width	16"-22"
Seat Depth	16"-22"
Backrest Height	21"-23" (18"-28" Proform)
Overall Weight	384 lbs (Power Tilt, Power AFP, 22NF Batts)
Controller	PG VR2-A2 60A
Warranty	18 Months: Motor and Electronics 5 year/ limited lifetime frame 1 Year Battery

HCPCS K0837

 Merits Health Products

 Meritshealthusa

 Merits USA

OPTIONS

- 22NF battery upgrade
- Power articulating foot platform
- Elevating leg rests or swing-away footrests
- Swing-away joystick
- USB charger



Merits Health Products Inc.

4245 Evans Avenue, Fort Myers FL, 33901

239-239-0579 | sales@meritsusa.com

merits