

Atlantis Specifications

Product Number	P710
Weight Capacity	600 lbs
Drivetrain	RWD
Seat Width	24 - 26"
Seat Depth	20 - 22"
Backrest Height	24"
Seat-to-Floor Height	23 - 24"
Overall Length	44"
Overall Width	27"
Overall Height	50"
Turning Radius	31.5"
Ground clearance	3"
Total Weight	304 lbs
Weight of Heaviest Piece	167 lbs
Max Speed	5 mph
Battery Size (2)	12V / Group 24 (75AH)
Per Charge Range	up to 32 miles
Motor	DC 24V, 300W
Controller	P&G VR2 60A
Charger	6A off-board
Gradient	10°
Front Wheels	9" - pneumatic
Drive Wheels	12" - pneumatic
Rear Wheels	3" - anti-tippers
HCPCS Code	K0826 / K0827
Standard Colors	Candy Apple Red or Metallic Blue

