

Axcel Narrow Specifications

Product Number	P335
Weight Capacity	300 lbs
Drivetrain	FWD
Seat Width	14-20"
Seat Depth	14-20"
Backrest Height	18-22"
Seat-to-Ground Height	17.5" - 19.5"
Seat- Elevation (Optional)	12"
Overall Length	39.6"
Overall Width	22.5"
Overall Height	41.5"
Turning Radius	26" without footplate
Ground clearance	3"
Total Weight	350-385 lbs
Max Speed	5.2 mph
Battery Size (2)	MK 50-12 SLD M 50AH
Per Charge Range	up to 19 miles
Motor	2-Pole / 222W
Controller	R-NET 90A
Charger	6 Amp
Gradient	7.5°
Front Wheels	Anti-tipper
Drive Wheels	14" x 3" hybrid gel
Rear Wheels	9" hybrid gel
HCPCS Code	K0856 / K0861
Standard Colors	Black Noir, Dove Gray, Sapphire Blue, Ruby Red, and Camo Green

