



# EZ-GO

K0814

*Enhancing life through mobility.*



Scan to learn more!

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features

- Ergonomic mid-back seat for additional back support
- Quick connect feature for easy assembly and transportation
- Turning radius of 23" for easy maneuverability in tight spaces
- Weight capacity, 300 lbs.
- Compact and lightweight design, weighs a total of 113 lbs.

## Accessories

- Oxygen tank holder
- Cup holder
- Cane & crutch holder
- Mounting bracket for cup holder
- Walker holder
- Side bag
- Foldable rear basket
- USB charger



## Specifications

Product#	P321
Weight Capacity	300 lbs.
Drivetrain	RWD
Seat Type	Stadium
Seat Width	18"
Seat Depth	18"
Backrest Height	17"
Seat-to-Floor Height	18" - 20"
Overall Length	27"
Overall Width	23"
Overall Height	43"
Turning Radius	23"
Ground Clearance	1.8"
Total Weight	113 lbs.
Wt of Heaviest Pc.	54 lbs.
Max Speed	4.5 mph
Range Per Charge	up to 12 mi.
Motor	DC 24V / 90W
Controller	P&G NVR2
Battery Size (2 required)	12V / 22AH
Charger	1.5A off-board
Gradient	6°
Front Wheels	7" - solid
Drive Wheels	8" - foam-filled
Rear Wheels	3" anti-tippers
HCPCSCode	K0814
Standard Colors	Red or Black

## Color Options



RED



BLACK



merits

Avid

Pilot