

## **Gemini Specifications**

---

Product Number	P3017 and P3011 with Lift
Weight Capacity	450 lbs and 350 lbs with Lift
Drivetrain	RWD
Seat Width	20 - 22"
Seat Depth	20 - 22"
Backrest Height	22"
Seat-to-Floor Height	22 - 29"
Seat Elevation	5" P3011 with Lift
Overall Length	41"
Overall Width	24"
Overall Height	45"
Turning Radius	21"
Ground clearance	4"
Total Weight	242 lbs and 350 lbs
Weight of Heaviest Piece	126 lbs and 146 lbs
Max Speed	5 mph
Battery Size (2)	12V / 22NF (50AH)
Per Charge Range	up to 20 miles
Motor	DC 24V, 240W
Controller	PG VR2 60A
Charger	5A off-board
Gradient	10°
Front Wheels	9" - foam-filled
Drive Wheels	14" - foam-filled
Rear Wheels	3" anti-tippers
HCPCS Code	K0824 / K0825 / K0830 / K0831
Standard Colors	Candy Apple Red or Metallic Silver Blue

