

## **Pioneer 10 Specifications**

---

Product Number	S341
Weight Capacity	500 lbs
Seat Width	22"
Seat Depth	20"
Backrest Height	22"
Seat-to-Ground Height	25" - 28"
Seat-to-Deck Height	16" - 19"
Overall Length	59"
Overall Width	27"
Overall Height	47"
Turning Radius	62"
Ground clearance	4"
Weight of Heaviest Piece	121 lbs
Total Weight	269 lbs
Max Speed	6 mph
Per Charge Range	up to 32 miles
Motor	DC 24V, 400W
Controller	Dynamic Rhino2 160A
Battery Size (2)	12V/Group 24 (70-80AH)
Charger	6A off-board
Gradient	12°
Front Wheel(s)	4-Wheel: 12.5" flat-free
Rear Wheels	12.5" flat-free
Suspension	Yes
HCPCS Code	K0899
Standard Color	Red

