

Pioneer 4 Specifications

Product Number	S141
Weight Capacity	400 lbs
Seat Width	20"
Seat Depth	20"
Backrest Height	20"
Seat-to-Ground Height	25" - 28"
Seat-to-Deck Height	16" - 19"
Overall Length	48"
Overall Width	23"
Overall Height	45.5"
Turning Radius	52"
Ground clearance	4"
Weight of Heaviest Piece	97 lbs
Total Weight	189 lbs
Max Speed	5 mph
Battery Size (2)	12V/U1(34AH)
Per Charge Range	up to 18 miles
Motor	DC 24V, 250W
Controller	Dynamic Rhino2 90A
Charger	5A off-board
Gradient	8°
Front Wheels	10" flat-free
Rear Wheels	10" flat-free
HCPCS Code	K0899
Standard Colors	Blue and Red

