

Pioneer 9 Specifications

Product Number	S331
Weight Capacity	500 lbs
Seat Width	22"
Seat Depth	20"
Backrest Height	22"
Seat-to-Ground Height	25" - 28"
Seat-to-Deck Height	16" - 19"
Overall Length	54"
Overall Width	27"
Overall Height	47"
Turning Radius	50"
Ground clearance	4"
Weight of Heaviest Piece	101 lbs
Total Weight	249 lbs
Max Speed	6 mph
Per Charge Range	up to 32 miles
Motor	DC 24V, 400W
Controller	Dynamic Rhino2 160A
Battery Size (2)	12V/Group 24 (70-80AH)
Charger	6A off-board
Gradient	12°
Front Wheel(s)	3-Wheel: 10" flat-free
Rear Wheels	12.5" flat-free
Suspension	Yes
HCPCS Code	K0808
Standard Colors	Red

