

Regal Specifications

Product Number	P310
Weight Capacity	300 lbs
Drivetrain	RWD
Seat Width	16, 18 and 20"
Seat Depth	16-20"
Backrest Height	20"
Seat-to-Floor Height	20-22"
Overall Length	39"
Overall Width	24"
Overall Height	41"
Turning Radius	21"
Ground clearance	3"
Total Weight	189 lbs
Weight of Heaviest Piece	97 lbs
Max Speed	5 mph
Battery Size (2)	12V / U1 (34AH)
Per Charge Range	up to 18 miles
Motor	DC 24V, 180W
Controller	P&G VR2 50A
Charger	5A off-board
Gradient	12°
Front Wheels	8" foam filled
Drive Wheels	10" foam filled
Rear Wheels	3" anti-tippers
HCPCS Code	K0822 / K0823
Standard Colors	Red or Blue

