



# VISION SUPER

K0824 / K0825 / K0830 / K0831

*Enhancing life through mobility*



Scan to learn more!

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features

- True mid-wheel: 20" turning radius
- Robust independent suspension can absorb up to 350 lbs. on the footplate without tipping forward (does not apply to lift option)
- Dynamic LiNX & Shark controller for convenient Bluetooth programming and diagnostics
- Premium, semi-reclining ergonomic seat up to 135°
- 20" - 22" customizable seat depth

## Accessories

- Unoccupied transit tie down kit
- Cup holder
- USB charger
- Oxygen tank holder
- Walker holder
- Foldable rear basket
- Safety flag
- Side under, arm bag

## Options

- Optional, 20" - 24" seat width
- 10" scissor power seat lift
- Elevating leg rests (K0195)
- Captain or Pan powered seat option

## Color Options



Red



Blue



## Specifications

Product#	P3277, P3274 w/lift
Weight Capacity	450 lbs., 400lbs. w/lift
Drivetrain	MWD
Seat Type	Captain or Pan
Seat Width	20" - 22"
Seat Depth	20" - 22"
Backrest Height	22"
Seat-to-Floor Height	22" - 34"
Seat Elevation	10" - P3274 w/lift
Overall Length	35.4"
Overall Width	27.2"
Overall Height	48"
Turning Radius	21"
Ground Clearance	2.2"
Total Weight	243 lbs., 273 lbs.
Wt of Heaviest Pc.	121 lbs., 151 lbs.
Max Speed	5 mph
Range Per Charge	up to 20 mi.
Motor	DC 24V / 300W x 2pcs
Controller	Dynamic Shark 50A / PG VR2 60A (Lift)
Battery Size (2 required)	12V / 22NF (50AH)
Charger	5A off-board
Gradient	8°
Front Wheels	6" - solid
Drive Wheels	12" - flat-free
Rear Wheels	6" - solid
HCPCS Codes	K0824 / K0825 K0830 / K0831
Standard Colors	Red or Blue

**merits**

merits *Pilot*

