

Roadster 4 Specifications

Product Number	S741A
Weight Capacity	300 lbs
Seat Width	18"
Seat Depth	17"
Backrest Height	15"
Seat-to-Ground Height	19" - 21"
Seat-to-Deck Height	16" - 18"
Overall Length	30"
Overall Width	21"
Overall Height	36"
Turning Radius	43.3"
Ground clearance	2"
Weight of Heaviest Piece	36 lbs
Total Weight	106 lbs
Max Speed	4 mph
Per Charge Range	up to 6 miles
Motor	DC 24V, 100W
Controller	PG S-DRIVE 45A
Battery Size (2)	12V/12AH
Charger	1.8A off-board
Gradient	6°
Front Wheel(s)	8" PU
Rear Wheels	8" PU
Suspension	n/a
HCPCS Code	K0800
Standard Color	Blue

