

## Travel-Ease HD Specifications

---

Product Number	P181
Weight Capacity	450 lbs
Drivetrain	RWD
Seat Width	22"
Seat Depth	18"
Backrest Height	15"
Seat-to-Floor Height	20"
Overall Length	45"
Overall Width	28"
Overall Height	36"
Turning Radius	25"
Ground clearance	4.8"
Total Weight	182 lbs
Weight of Heaviest Piece	97 lbs
Max Speed	4 mph
Battery Size (2)	12V / 22NF (50AH)
Per Charge Range	up to 25 miles
Motor	DC 24V, 300W
Controller	VR2 70A
Charger	5A off-board
Gradient	6-10°
Front Wheels	9" - foam-filled
Drive Wheels	12.5" - foam-filled
Rear Wheels	1.5" anti-tippers
HCPCS Code	K0824
Standard Colors	Silvertone

