

Travel Ease Specifications

Product Number	P101
Weight Capacity	300 lbs
Drivetrain	RWD
Seat Width	16, 18 and 20"
Seat Depth	16"
Backrest Height	15-19"
Seat-to-Ground Height	20"
Overall Length	43"
Overall Width	22, 24, and 26"
Overall Height	36"
Turning Radius	23"
Ground clearance	4.8"
Total Weight	146, 150, and 154 lbs
Weight of Heaviest Piece	86 - 97 lbs
Max Speed	4 mph
Battery Size (2)	12V / U1 (34AH)
Per Charge Range	up to 18 miles
Motor	DC 24V, 180W
Controller	VSI
Charger	5A off-board
Gradient	10°
Front Wheels	8" foam-filled
Drive Wheels	12.5" foam-filled
Rear Wheels	1.5" anti-tippers
HCPCS Code	K0822
Standard Colors	Silvertone

