

Ultra CST Specifications

Product Number	P325
Weight Capacity	300 lbs
Drivetrain	MWD
Seat Width	18" or 20"
Seat Depth	18" or 20"
Backrest Height	20"
Seat-to-Ground Height	21.5" - 22.5.5"
Tilt (Optional)	0° - 50°
Recline	90° - 105° (manual)
Overall Length	36"
Overall Width	24"
Overall Height	40"
Turning Radius	28" including footplate
Ground clearance	2.5"
Total Weight	310 lbs
Max Speed	5 mph
Battery Size (2)	24V / 22NF
Per Charge Range	up to 11 miles / 19 miles
Motor	DC 24V, 240W
Controller	PG VR2-A2 60A
Charger	6 Amp
Gradient	6°
Front Wheels	6" x 2"
Drive Wheels	12" x 3"
Rear Wheels	6" x 2"
HCPCS Code	K0836
Standard Colors	Red

