



# ULTRA

K0822 / K0823 / K0835

*Enhancing life through mobility.*



Scan to learn more!

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

©Merits Health Products, Inc. rev. 1024

## Features

- Center of gravity tilt 0° - 50° posterior
- Bright LED lights, included
- Shear reduction recline 90° - 150°
- Tension-adjustable front and rear suspension
- Adjustable seat width and depth from 16" - 22"; depth adjustable in 1" increments
- 24" base width and seat-to-floor height as low as 16.5"
- Rated up to 18 miles per charge

## Options

- 22NF battery upgrade
- 12" elevating seat
- Power articulating foot platform
- Elevating leg rests or swing-away footrests
- Swing-away joystick
- Oxygen tank holder
- USB charger

## Models

Vision Ultra - G2, Standard, Captain P325-2S C

Vision Ultra - G2, Standard, Pan P325- 2S S

Vision Ultra - G2, Standard, Single Power P325-2SSP S



## Specifications

Product#	P325
Weight Capacity	300 lbs.
Drivetrain	MWD
Maximum Speed	5 mph
Per-Charge Range	up to 19 mi.
Battery Size (2 required)	35Ah, U1; 55Ah, 22NF, 55Ah, 22NF GEL
Ground Clearance	2.5"
Turning Radius	26"
Overall Length	36"
Overall Width	24"
Suspension	Front and Rear
Seat-to-Floor Height	Power tilt or pan: 16.5"-19.5", Captain: 20"-23"
Seat Width	16"-22", adjusted in 1" increments
Seat Depth	16"-22", adjusted in 1" increments
Backrest Height	18"-28"
Total Weight	315 - 384 lbs.
Motor	2-Pole/240W
Controller	P&G VR2-A2 60A
Charger	6 Amp
Gradient	6°
Front Wheels	Anti-Tipper
Drive Wheels	14"x 3" foam-filled
Casters	6" x 2" foam-filled

## Color Options



White

Red

Blue

Charcoal  
Gray

Optional MK batteries



**+** **merits**

merits *Avid* *Pilot*

