

Ultra HD Specifications

Product Number	P331
Weight Capacity	450 lbs
Drivetrain	MWD
Seat Width	20-24"
Seat Depth	16-20"
Backrest Height	21-23"
Seat-to-Ground Height	17.5" - 18.5"
Tilt (Optional)	0° - 50°
Recline (Optional)	90° - 110° fixed
Articulating Footplate	Optional
Overall Length	35" / 45" with footplate
Overall Width	27"
Overall Height	42"
Turning Radius	21"
Ground clearance	2.875"
Total Weight	288-310 lbs
Max Speed	5 mph
Battery Size (2)	34Ah/U1 or 55Ah/22NF
Per Charge Range	up to 8 miles / 19 miles
Motor	2-Pole/200W
Controller	PG VR2-A2 60A
Charger	6 Amp
Gradient	6°
Front Wheels	6" x 2"
Drive Wheels	12" x 3"
Rear Wheels	6" x 2"
HCPCS Code	K0837
Standard Colors	Red, Blue, White, and Gray

