

## Ultra Specifications

---

Product Number	P325
Weight Capacity	300 lbs
Drivetrain	MWD
Seat Width	16-22"
Seat Depth	16-22"
Backrest Height	18-28"
Seat-to-Ground Height	16.5" - 19.5"
Seat- Elevation (Optional)	12"
Tilt (Optional)	0°-50°
Recline (Optional)	90°-110° fixed
Articulating Footplate	Optional
Overall Length	36"
Overall Width	24"
Overall Height	40"
Turning Radius	26"
Ground clearance	2.5"
Total Weight	315-384 lbs
Max Speed	5 mph
Battery Size (2)	34Ah/U1 or 55Ah/22NF
Per Charge Range	up to 12 miles / 19 miles
Motor	2-Pole/240W
Controller	PG VR2-A2 60A
Charger	6 Amp
Gradient	6°
Front Wheels	6" x 2"
Drive Wheels	14" x 3"
Rear Wheels	6" x 2"
HCPCS Code	K0822 / K0823 / K0835
Standard Colors	Red, Blue, Charcoal Gray and White

