

Ultra X Specifications

| | |
|---------------------------|----------------------------|
| Product Number | P325 |
| Weight Capacity | 330 lbs |
| Drivetrain | MWD |
| Seat Width | 16-22" |
| Seat Depth | 16-22" |
| Backrest Height | 18-28" |
| Seat-to-Ground Height | 16.5-19.5" |
| Elevating Seat (Optional) | 12" |
| Tilt (Optional) | 0° - 50° |
| Recline (Optional) | 90° - 110° fixed |
| Articulation Footplate | Optional |
| Overall Length | 36" |
| Overall Width | 24" |
| Overall Height | 40" |
| Turning Radius | 26" |
| Ground clearance | 2.5" |
| Total Weight | 315-384 lbs |
| Max Speed | 5 mph |
| Battery Size (2) | 34Ah/U1 or 55Ah/22NF |
| Per Charge Range | up to 12 miles / 19 miles |
| Motor | 2-Pole/240W |
| Controller | P&G VR2-A2 60A |
| Charger | 6 Amp |
| Gradient | 6° |
| Front Wheels | 6 x 2" |
| Drive Wheels | 14 x 3" |
| Rear Wheels | 6 x 2" |
| HCPCS Code | K0837 |
| Standard Colors | Red, Blue, White, and Gray |

