

Vision Sport Specifications

Product Number	P326E
Weight Capacity	300 lbs
Drivetrain	MWD
Seat Width	16, 18 and 20"
Seat Depth	16 - 20"
Backrest Height	20"
Seat-to-Floor Height	17.5, 19.5 and 21.5"
Overall Length	34"
Overall Width	24"
Overall Height	48.5"
Turning Radius	20"
Ground clearance	2"
Total Weight	184 lbs
Weight of Heaviest Piece	90 lbs
Max Speed	5 mph
Battery Size (2)	12V / U1 (34AH)
Per Charge Range	up to 18 miles
Motor	DC 24V, 160W
Controller	Dynamic LiNX 40A
Charger	4A off-board
Gradient	6°
Front Wheels	6" solid
Drive Wheels	10" solid
Rear Wheels	6" solid
HCPCS Code	K0822 / K0823
Standard Colors	Red and Blue

