

Vision Super Specifications

Product Number	P3277 and P3274 with Lift
Weight Capacity	450 lbs and 400 lbs with Lift
Drivetrain	MWD
Seat Width	20- 22"
Seat Depth	20 - 22"
Backrest Height	22"
Seat-to-Floor Height	22 - 34"
Seat Elevation	10" P3274 with Lift
Overall Length	35.4"
Overall Width	27.2"
Overall Height	48"
Turning Radius	21"
Ground clearance	2.2"
Total Weight	243 lbs and 273 lbs
Weight of Heaviest Piece	121 lbs and 151 lbs
Max Speed	5 mph
Battery Size (2)	12V / 22NF (50AH)
Per Charge Range	up to 20 miles
Motor	DC 24V, 300W x 2 pieces
Controller	Dynamic Shark 50A / PG VR2 60A with Lift
Charger	5A off-board
Gradient	8°
Front Wheels	6" - solid
Drive Wheels	12" - flat-free
Rear Wheels	6" - solid
HCPCS Code	K0824 / K0825 / K0830 / K0831
Standard Colors	Red or Blue

