




GEMINI

K0824 / K0825 / K0830 | K0831

Enhancing life through mobility



Scan to learn more!

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- Ergonomic semi-reclining seat that reclines up to 135°
- 14" rubber drive tires for superior traction
- Dynamic LiNX controller for easy programming and diagnostics
- 4" ground clearance for easy obstacle clearance
- Optional 5" power elevating seat.

Accessories

- Oxygen tank holder
- Cane & crutch holder
- Walker holder
- Foldable rear basket
- Cup holder
- Side bag
- Cover
- USB charger

Options

- Captain or Pan powered seat option
- 20" and 22" seat width option
- Elevating leg rest available (K0195)
- 5" power lift option, available
- Swing-away joystick
- XL footplate

Color Options



Candy Apple Red



Metallic Silver Blue



Specifications

Product#	P3017, P3011 w/lift
Weight Capacity	450 lbs., 350 lbs. w/lift.
Drivetrain	RWD
SeatType	Captain or Pan
Seat Width	20" - 22"
Seat Depth	20" - 22"
Backrest Height	22"
Seat-to-Floor Height	22" - 29"
Seat Elevation	5" - P3011
Overall Length	41"
Overall Width	24"
Overall Height	45"
Turning Radius	21"
Ground Clearance	4"
Total Weight	242 lbs., 350 lbs.
Wt of Heaviest Pc.	126 lbs., 146 lbs.
Max Speed	5 mph
Range Per Charge	up to 20 mi.
Motor	DC 24V, / 240W
Controller	P&G VR2 60A
Battery Size (x2)	12V / 22NF (50AH)
Charger	5A off-board
Gradient	10°
Front Wheels	9" - foam-filled
Drive Wheels	14" - foam-filled
Rear Wheels	3" anti-tippers
HCPCS codes	K0824 / K0825 K0830 / K0831

