




VECTOR HD

K0850 / K0858 / K0862

Enhancing life through mobility.



Contact your Territory Sales Manager or email us at sales@meritsusa.com

Scan to learn more!

©Merits Health Products, Inc. rev. 1024

Features

- Greater stability with wider casters and weight capacity of up to 450 lbs.
- Pre-drilled for easy installation, saving you hundreds of dollars
- Full front and rear suspension to tackle all those bumps and curbs
- Dual in-line motors for enhanced efficiency, torque, range and performance
- Rated up to 12.5 miles per charge

Accessories

(Based on chair configuration)

- USB charger
- Oxygen tank holder

Options

- Joystick knob choices
- Multiple armrest add-ons
- Multiple footplate options available
- Joystick mounting
- Captain, Pan or AVID powered seat



Optional MK batteries



Specifications

Product#	P323
Weight Capacity	450 lbs.
Drivetrain	MWD
Seat Width	20" - 24"
Seat Depth	16" - 22"
Backrest Height	18" - 28"
Seat-to-Floor Height	18.5" - 19.5"
Seat Elevation (Optional)	12"
Overall Length	47"
Overall Width	25"
Overall Height	42"
Turning Radius	27"
Ground Clearance	3"
Total Weight	315 - 416 lbs.
Max Speed	5 mph
Battery Size (2 required)	55 AH, 22 NF Gel, 22 MF, AGM
Range Per Charge	up to 12.5 mi.
Motor	4-Pole / 200W
Controller	P&G VR2 or R-NET
Charger	6 Amp
Gradient	6°
Front Wheels	8" x 3" foam-filled
Drive Wheels	14" x 3" foam-filled
Rear Wheels	8" x 3" foam-filled
HCPCS Code	K0850 / K0858 K0862

Color Options

