




VISION SPORT

K0822 / K0823

Enhancing life through mobility



Scan to learn more!

Contact your Territory Sales Manager or email us at sales@meritsusa.com

Features

- True mid-wheel: 20" turning Radius
- Robust independent suspension can absorb up to 250 lbs. on the footplate without tipping.
- Premium, semi-reclining ergonomic seat, up to 135°
- 10" Solid drive tires for superior traction
- Dynamic LiNX 40A controller for convenient Bluetooth programming and diagnostics
- 16" - 20" customizable seat depth

Options

- Unoccupied transit tie down kit
- Cup holder
- USB charger
- Oxygen tank holder
- Walker holder
- 16" - 20" seat width option
- Swing-away joystick
- Elevating leg rests (K0195)
- Foldable rear basket
- Cane & crutch holder



Specifications

Product#	P326
Weight Capacity	300 lbs.
Drivetrain	MWD
Seat Type	Captain or Pan
Seat Width	16" - 20"
Seat Depth	16" - 20"
Backrest Height	18" - 20"
Seat-to-Floor Height	17.5" - 23"
Overall Length	34"
Overall Width	24"
Overall Height	48.5"
Turning Radius	20.6"
Ground Clearance	2"
Total Weight	184 lbs.
Wt of Heaviest Pc.	90 lbs.
Max Speed	5 mph
Range Per Charge	up to 18 mi.
Battery Size (2 required)	12V/U1 (34AH)
Motor	DC24V / 160W
Controller	Dynamic LiNX 40A
Charger	4A off-board
Gradient	6°
Front Wheels	6" - solid
Drive Wheels	10" - solid
Rear Wheels	6" - solid
HCPCSCode	K0822 / K0823
Standard Colors	Red or Blue

Color Options



Red



Blue

