

merits

# VISION SPORT P326D

MWD Powerbase with Lift

K0830 / K0831



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

# Features

- True mid-wheel: 20" turning radius
- Toggle lift operation or upgrade through the joystick (E2310)
- VR2 joystick is programmed via PP1 programmer
- Premium, semi-reclining ergonomic seat up to 135°
- 10" solid tires for superior traction
- 16" - 20" customizable seat depth
- 7.5" power elevating seat



# Accessories

- Unoccupied transit tie down kit
- Cup holder
- USB charger
- Oxygen tank holder
- Walker holder
- Foldable rear basket
- Safety flag
- Side bag

# Options

- Optional, 16" - 20" seat width
- Swing-away joystick
- Elevating leg rests (K0195)
- Captain or Pan powered seat option

# Specifications

Product#	P326D
Weight Capacity	300 lbs.
Drive train	MWD
Seat Type	Captain or Pan
Seat Width	18" - 20"
Seat Depth	16" - 20"
Backrest Height	20"
Seat-to-Floor Height	20.5" - 29.5"
Seat Elevation	7.5"
Overall Length	34"
Overall Width	24"
Overall Height	48.5"
Turning Radius	20"
Ground Clearance	2"
Total Weight	214 lbs.
Wt of Heaviest Pc.	120 lbs.
Max Speed	5 mph
Range Per Charge	up to 18 mi.
Motor	DC 24V / 160W
Controller	P&G VR2 60A
Battery Size (2 required)	12V / U1 (34AH)
Charger	4A off-board
Gradient	6°
Front Wheels	6" - solid
Drive Wheels	10" - solid
Rear Wheels	6" - solid
HCPCS Code	K0830 / K0831
Color Options	 Red  Black



 **Enhancing life through mobility.**

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com

