

Dualer Specifications

Product Number	P312
Weight Capacity	300 lbs
Drivetrain	RWD/FWD
Seat Type	Captain or Pan
Seat Width	16" - 18"
Seat Depth	18"
Backrest Height	20"
Seat-to-Floor Height	21.5" - 26"
Elevating Seat	5"
Overall Length	38"
Overall Width	24"
Overall Height	41"
Turning Radius	21"
Ground clearance	2"
Total Weight	189 lbs
Weight of Heaviest Piece	97 lbs
Max Speed	5 mph
Battery Size (2)	12V / U1 (34AH)
Per Charge Range	up to 18 miles
Motor	DC 24V / 180W
Controller	P&G VR2 60A
Charger	5A off-board
Gradient	12°
Front Wheels	8" - solid
Drive Wheels	10" - foam-filled
Rear Wheels	3" - anti-tippers
HCPCS Code	K0830 / K0831
Standard Colors	Candy Apple Red

