

merits

**DUALER** P312

RWD / FWD Powerbase

K0830 / K0831



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

# Features

- The Dualer can change from RWD to FWD with swivel of the seat
- Ergonomic semi-reclining seats that recline up to 135°
- 10" knobby drive tires for superior traction
- Dual footplates for comfort and convenience
- 5" power elevating seat.

# Accessories

- Oxygen tank holder
- Cane & crutch holder
- Walker holder
- Foldable rear basket
- Cup holder
- Side bag
- Cover

# Options

- Captain or Pan powered seat option
- Swing-away joystick
- 16" and 18" seat width option
- Elevating legrest available (K1095)

# Specifications

Product#	P312
Weight Capacity	300 lbs.
Drivetrain	RWD / FWD
SeatType	Captain or Pan
Seat Width	16" or 18"
Seat Depth	18"
Backrest Height	20"
Seat-to-Floor Height	21.5" - 26"
Seat Elevation	5"
Overall Length	38"
Overall Width	24"
Overall Height	41"
Turning Radius	21"
Ground Clearance	2"
Total Weight	189 lbs.
Wt of Heaviest Pc.	97 lbs.
Max Speed	5 mph
Range Per Charge	up to 18 mi.
Motor	DC24V / 180W
Controller	P&G VR2 60A
Battery Size (x2)	12V / U1 (34AH)
Charger	5A off-board
Gradient	12°
Front Wheels	8" - solid
Drive Wheels	10" - foam-filled
Rear Wheels	3" anti-tippers
HCPCS Codes	K0830 / K0831
Color Options	<span style="color: red;">■</span> Candy Apple Red



merits



**merits** *Enhancing life through mobility.*

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0225