

Ez-Go Specifications

Product Number	P321
Weight Capacity	300 lbs
Drivetrain	RWD
Seat Type	Stadium
Seat Width	18"
Seat Depth	18"
Backrest Height	17"
Seat-to-Ground Height	18-20"
Overall Length	27"
Overall Width	23"
Overall Height	43"
Turning Radius	23"
Ground clearance	1.8"
Total Weight	113 lbs
Weight of Heaviest Piece	54 lbs
Max Speed	4.5 mph
Battery Size (2)	12V / 22AH
Per Charge Range	up to 12 miles
Motor	DC 24V / 90W
Controller	P&G NVR2
Charger	1.5A off-board
Gradient	6°
Front Wheels	7" solid
Drive Wheels	8" foam-filled
Rear Wheels	3" anti-tippers
HCPCS Code	K0814
Standard Colors	Red or Black

