

Pioneer 3 Specifications

Product Number	S131
Weight Capacity	400 lbs
Drivetrain	RWD
Seat Width	20"
Seat Depth	20"
Backrest Height	20"
Seat-to-Ground Height	25" - 28"
Seat-to-Deck Height	16" - 19"
Overall Length	47"
Overall Width	23"
Overall Height	45.5"
Turning Radius	39"
Ground clearance	4"
Weight of Heaviest Piece	86 lbs
Total Weight	178 lbs
Max Speed	5 mph
Per Charge Range	up to 18 miles
Motor	DC 24V, 250W
Controller	Dynamic Rhino2 90A
Battery Size (2)	12V/U1 (34AH)
Charger	5A off-board
Gradient	8°
Front Wheel(s)	10" flat-free
Rear Wheels	10" flat-free
Suspension	Yes
HCPCS Code	K0807
Standard Colors	Red

