

merits

# PIONEER 3

S131

3 Wheel Scooter

K0807



***Enhancing Life Through Mobility.***

Contact your Territory Sales Manager or email us at [sales@meritsusa.com](mailto:sales@meritsusa.com)

## Features

- Hydraulic, adjustable tiller
- 20"W x 20"D x 20"H, sliding, captain seat
- Half speed inhibit switch
- Headlight, tail lights, turn signals and hazard lights
- 4" ground clearance
- Three wheels provide a tighter turning radius for increased maneuverability

## Accessories

- Oxygen Tank Holder
- Cane & Crutch Holder
- Walker Holder
- (Lg) Foldable Rear Basket
- (Med) Foldable Rear Basket
- Cover

## Specifications

Product#	S131
Weight Capacity	400 lbs.
Drivetrain	RWD
Maximum Speed	5 mph
Per-Charge Range	up to 18 miles
Battery Size (2 required)	12V/U1 (34AH)
Ground Clearance	4"
Turning Radius	39"
Overall Length	47"
Overall Width	23"
Overall Height	45.5"
Suspension	Yes
Seat-to-Floor Height	16" - 19"
Seat-to-Ground Height	25" - 28"
Seat Width	20"
Seat Depth	20"
Backrest Height	20"
Total Weight	178 lbs.
Weight of Heaviest Pc.	86 lbs.
Motor	DC 24V, 250W
Controller	Dynamic Rhino2 90A
Charger	5A off-board
Gradient	8°
Front Wheels	10" flat-free
Drive Wheels	10" flat-free
HCPCS Code	K0807
Color Option	<span style="color: red;">■</span> Red



 **merits** *Enhancing life through mobility.*

©Merits Health Products, Inc.  
4245 Evans Avenue Fort Myers, FL 33901  
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



rev. 0225

merits