

merits

PIONEER 4 S141

Four Wheel Scooter w/Lift Option

K0899



Enhancing Life Through Mobility.

Contact your Territory Sales Manager or email us at sales@meritsusa.com



Features

- One-piece frame holds up to 400 lbs.
- Steel bumper comes standard
- Adjustable captain seat
- Padded foam armrests are height, weight, and angle adjustable
- Hazard lights
- Turn signals
- Headlight and tail lights
- Anti-tippers
- Range of up to 18 miles per charge
- 4" ground clearance
- Four wheels provide a tighter turning radius for increased maneuverability
- Optional seat lift

Accessories

- Oxygen tank holder - attaches to seat bracket
- Cane & crutch holder - attaches to seat bracket
- Walker holder - attaches to seat bracket
- Foldable rear basket (medium or large) - attaches to seat bracket

Specifications

Product#	S141
Weight Capacity	400 lbs
Drivetrain	RWD
Maximum Speed	5 mph
Per-Charge Range	up to 18 miles
Battery Size (2 required)	12V/U1(34AH)
Ground Clearance	4"
Turning Radius	52"
Overall Length	48"
Overall Width	23"
Overall Height	45.5"
Seat-to-Ground Height	25" - 28"
Seat-to-Deck Height	16" - 19"
Seat Width	20"
Seat Depth	20"
Backrest Height	20"
Total Weight	189 lbs
Motor	DC 24V, 250W
Controller	Dynamic Rhino2 90A
Charger	5A off-board
Gradient	8°
Front Wheels	10" flat-free
Drive Wheels	10" flat-free
HCPCS Code	K0899
Color Options	 Red  Blue



 **Enhancing life through mobility.**

©Merits Health Products, Inc.
4245 Evans Avenue Fort Myers, FL 33901
www.meritsusa.com | 239-772-0579 | sales@meritsusa.com



merits

rev. 0225