

Roadster S3 Specifications

Product Number	S731A
Weight Capacity	300 lbs
Drivetrain	RWD
Seat Width	18"
Seat Depth	17"
Backrest Height	15"
Seat-to-Ground Height	19" - 21"
Seat-to-Deck Height	16" - 18"
Overall Length	40"
Overall Width	21"
Overall Height	36"
Turning Radius	35.4"
Ground clearance	2"
Weight of Heaviest Piece	29 lbs
Total Weight	99 lbs
Max Speed	4 mph
Per Charge Range	up to 6 miles
Motor	DC 24V, 100W
Controller	PG S-DRIVE 45A
Battery Size (2)	12V/12AH
Charger	1.8A off-board
Gradient	6°
Front Wheel(s)	7" PU
Rear Wheels	8" PU
Suspension	n/a
HCPCS Code	K0800
Standard Colors:	Red and Blue

